




SamanTree  
Medical SA

# Histolog<sup>®</sup> Scanner

Cancer cells  
at your fingertips

Prostate brochure 2022



Robot-assisted radical prostatectomy (RARP) allows surgeons to remove the prostate with precision, making conservative treatment possible. Nerve-bundles preservation is now often offered to patients to maximize their chance of remaining continent for urine and maintaining erectile function.



# The Clinical Need

The absence of adequate intra-operative margin assessment (IOA) makes nerve-sparing surgeries prone to recurrences, with up to 40% of patients detected with positive margins<sup>1</sup>.

Neurovascular structure-adjacent frozen-section examination (NeuroSAFE) was introduced in 2005 to enable assessment of posterolateral margins adjacent to the neurovascular bundles. This technique enables a precise IOA, however, it is an elaborated procedure requiring trained resources and taking up to 1 hour. Only a few centers have set up this procedure as their clinical routine.

Prerequisites to standardization of nerve-sparing approach :

- Accurate IOA
- Cost-effective technique
- Technology usable by both surgeons and pathologist


<sup>1</sup>Jeong et al., 2022, Izard et al., 2015, Iczkowski et al., 2011



# Our mission

SamanTree Medical aims to be a game changer in the era of clinical workflow digitalization. We are committed to improving the journey of patients suffering from cancer by enabling fresh tissue analysis in real time to drastically reduce delays in establishing and executing the treatment plan.

---



### ENSURE SUPERIOR PATIENT CARE

Nerve-sparing approach is associated with an increased risk that tumour will be left on the surface of the resected specimen (positive margins).

### SAVE PRECIOUS TIME AND RESOURCES

Frozen section analysis is performed to evaluate the nerve-bundles. It requires at least a technician and a pathologist, resulting in a approx. 1 hour procedure.

### STANDARDIZATION OF NERVE-SPARING APPROACH

The risk of recurrence and IOA complex logistics make nerve-sparing available for only a subset of patients.

PROBLEM	High rate of positive margins associated to nerve-sparing approach	Frozen section analysis is time/ resource consuming	Nerve-sparing associated to IOA is proposed to only a subset of patients
SOLUTION	Real-time morphology information of the whole prostate, including areas adjacent to the neurovascular bundles.	Easy-to-use digital microscope for fresh tissue imaging in less than a minute. Minimal ressources utilization (1 clinician), compatible with current workflow (not damaging)	Providing a cost-effective solution to standardize the practice

# Our solution The Histolog® Scanner

## Global mapping of the prostate immediately during surgery

The Histolog® Scanner is a breakthrough medical imaging modality based on a novel ultra-fast confocal microscopy technology invented in 2010.

Its innovative design makes it highly practical for quick assessment during surgery, bringing the clinician one touch-on-the-screen away from visualizing cancerous cells immediately on a surgical specimen.

In RARP, an initial feasibility study showed that the Histolog Scanner enabled NeuroSAFE-like procedure with similar outcomes, however with > 80% time reduction and used by a pathologist alone.



A cost-efficient technique that provides real time morphology information



Easy-to-use platform: plug and play device with quick learning, usable by clinicians or OR staff



Digital images, enabling digital and remote workflows



Result in minutes : 15s for specimen preparation (10s fluorescent dye + 2s rinsing in saline solution) and ~50s for full-resolution image (large field of view 4.8cm x 3.6cm)



A dedicated team and network to support you in the implementation



# QUICK & CLEAN

4-steps procedure for  
accurate margin assessment  
immediately in the OR




## 1 Excision

Excise the tumor from the patient.



## 2 Preparation

Immerse the excision in Histolog Dip and rinse it.




## 4 Evaluation

With this additional information, return to patient in confidence.

Excision is sent for standard post-surgical pathology assessment.

### Touchscreen interface

Instant access to special features such as reporting & annotating tools




## 3 Imaging

Map in minutes the whole excision surface.

Excision remains visible and accessible during the entire imaging procedure.

**Histolog Dish**  
Single use receptacle



# Our partners

---

## Clinical partners

The Histolog Scanner is a game changing approach under evaluation by leading centers in Europe for its application in breast, pathology and prostate.

An initial feasibility study was performed with Canisius Wilhelmina Hospital (Netherlands), showing that the Histolog Scanner had similar performances as NeuroSAFE, however with drastic time savings (80%).

These promising data raised the possibility to having a cost-efficient technique for IOA in nerve-sparing surgery. The Histolog® Scanner is currently under evaluation in the course of the NeuroSAFE PROOF clinical trial initiated by the medical team of University College London (UCL) led by Mr Greg Shaw.

## Reference centers

Canisius Wilhelmina  
Hospital, Netherlands

University College London,  
United-Kingdom



King's College Hospital

Imperial College  
Hospital

**University College  
London Hospital**

Erasmus Hospital

**Canisius Wilhelmina Hospital**

St. Vincenz Hospital

Heidelberg University Hospital

University Hospital rechts der Isar

Brust Centrum Zürich

Gustave Roussy Institute

Institut Cancer de Lorraine

Valais Hospital

Institut Bergonié

Reggio Emilia Hospital


# Learn to use Histolog images in no-time

## Histolog Image Training program (HIT)

**Short. Flexible. Simple.**

**A learning program designed for clinicians.**

The Histolog Image Training (HIT) was developed with our community of pathologists and experts to provide a simple and efficient way of getting familiar with Histolog image. Designed for both beginners and experienced morphology content readers, the HIT is accessible to all and allows for flexible learning. In and out of the operating theatre, you set the pace and we keep it.



# THE GOLD STANDARD IN THE OR

Once resected, the full excision is simply dipped in a contrast agent and placed on a disposable dish over the imaging window. Within 1 minute, the specimen is ready for imaging. There is no need for fixing, nor embedding, nor glass mounting of the excision. Then, within a few seconds, a digital image of the excision surface is displayed on a touch screen which allows to zoom and move around in the image. Surgeon can visualize suspicious lesion at a glance. Importantly, the excision remains intact for later final assessment.



# Examples of peripheral area images close to neurovascular bundles

5% Zoom level (Full Field of View) of Histolog Scanner	<div data-bbox="89 129 740 163" data-label="Caption"><p>Tissue Sections with Vessel, Muscle and Normal Glands</p></div> <div data-bbox="185 203 632 768" data-label="Image"></div> <div data-bbox="97 781 772 873" data-label="Caption"><p>Large normal glands can be identified at this magnification. A large vessel is seen within a green annotation. The black frame indicates the area that is shown at the highest magnification on the right.</p></div>	<div data-bbox="850 129 1227 163" data-label="Caption"><p>Prostate cancer at the periphery</p></div> <div data-bbox="857 208 1543 689" data-label="Image"></div> <div data-bbox="861 781 1536 873" data-label="Caption"><p>Large normal glands can be identified at this magnification. A large vessel is seen within a green annotation. The black frame indicates the area that is shown at the highest magnification on the right.</p></div>
25% Zoom level (Full Field of View) of Histolog Scanner	<div data-bbox="97 931 780 1400" data-label="Image"></div> <div data-bbox="97 1426 772 1518" data-label="Caption"><p>Peripheral area close to neurovascular bundles. Large normal glands on the left, large vessel in the center and muscular fibers with connective tissue on the right.</p></div>	<div data-bbox="863 936 1549 1364" data-label="Image"></div> <div data-bbox="863 1393 1522 1518" data-label="Caption"><p>Small suspicious glandular structures are seen on the border of the specimen (artefact of organ dissection exposing inner structures on the surface of the specimen that should not be considered as positive margin implying intraoperative actions).</p></div>
100% Zoom level (Full Field of View) of Histolog Scanner	<div data-bbox="97 1588 780 2049" data-label="Image"></div> <div data-bbox="97 2072 424 2163" data-label="Caption"><p>Large normal glands on the left (<b>G</b>) Large vessel in the center (<b>V</b>) and Muscular fibers (<b>M</b>).</p></div>	<div data-bbox="863 1588 1549 2011" data-label="Image"></div> <div data-bbox="861 2072 1511 2163" data-label="Caption"><p>Cancerous glands are seen in the image. They are small roundish histological structures presenting an epithelial texture surrounded by stroma.</p></div>



Distributed by:



 [www.axlab.dk](http://www.axlab.dk)

 [axlab@axlab.dk](mailto:axlab@axlab.dk)

 +45 35 43 00 73



SamanTree  
Medical SA

© OCTOBER 2022 SamanTree Medical SA

All rights reserved