



## BD Totalys™ SlidePrep

Fully automated slide processing  
from preparation to staining\*

- Maintain **continuous chain of custody** through integration with the **BD Totalys™ MultiProcessor** and **DataLink middleware connection**
- **Random sample loading**
- **Automatic barcode scanning** ensures sample match and chain of custody
- Process **BD SurePath™ Liquid-based Pap** or **non-GYN samples**

A new level of confidence with  
automated individual staining

- **Consistent results** with staining in individual chambers using fresh reagents for every sample
- **Accurate** reagent pipetting and **precise** timing for each sample reduces process variability

\*When used with BD Totalys™ MultiProcessor





# Help improve your lab's productivity

## Efficient throughput

- Complete a single run of **48 slides** in **68 minutes**
- High throughput capability: **Up to 336 stained slides** produced from cell pellet **in an 8-hour shift**
- **Minimal hands-on-time** allows staff to perform other tasks

### Specimen handling

Specimen collection type	BD SurePath™ Collection Vial and CytoRich™ non-gynaecological samples
On-board sample capacity	BD SurePath™ Liquid-based Pap Test and CytoRich™ non-gynaecological tests: Up to 48 samples
Positive sample identification	2D barcodes on all centrifugation tubes and cytology slides ensure chain of custody throughout processing

### Dimensions

Height	86.4 cm (34 in)
Width (including bulk reagent station and disposable tip waste container)	182.9 cm (72 in)
Depth	62.2 cm (24.5 in)
Weight	124.7 kg (275 lbs)
Separate bulk reagent station	38 cm H x 35.5 cm W x 53.5 cm D (15" H x 14" W x 21" D)
Separate waste station dimensions	40.6 cm H x 33 cm W x 58.4 cm D (16" H x 13" W x 23" D)

### Operating conditions

Temperature	18 – 30 °C (64 – 86 °F)
Humidity	20 – 80% RH, non condensing
Altitude	Up to 2,000 m (6,562 ft)

### Electrical requirements

Input voltage	100 – 240 VAC
Rated current	12 A
Input line frequency	50 – 60 Hz
Power	1,000 Watts



BD Totalys™ SlidePrep  
Automated slide preparation and staining

Ready to see **a real difference?**

Contact your BD Life Sciences Account Executive today.

BD - Europe, Terre-Bonne Park - A4, Route de Crassier 17, 1262 Eysins, Switzerland

**bd.com**

Products CE Marked in compliance with the European In Vitro Diagnostic Medical Device Directive 98/79/EC.  
© 2018 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.  
XEUR5332-17

