



205300 Histo Bath HIR-3D

## Kunz Histo Bath HIR-3D

### Advantages

- LED-light gives excellent visibility in the water
- Large round warm area, drying up to 27 slides
- No need for extra heating plate
- Easy to use digital touch control
- Light weight
- Easy and convenient to carry and empty the water
- An optional lid is available to cover the water
- Kunz well-known product quality

### Specifications

Mains power:	220-240VAC 50-60Hz
Power consumption:	190W
Temperature range:	+20 to +80°C
Temp. accuracy:	± 1°C
Internal dimensions:	Diameter 217mm, Depth 61mm
External dimensions:	Diameter 323mm, Height 88mm
Weight:	1,5kg

