


# C-Scan

Automatic Cassette Detecting System




**MILESTONE**  
H E L P I N G  
P A T I E N T S

## AUTOMATED BIOPSIES CAPTURING



### | DIGITAL DOCUMENTATION OF BIOPSIES

Pathology laboratories deal with hundreds, often thousands, of biopsies every day. Documenting their preparation is a time-consuming task and requires adequate tools which are normally not integrated into the process. The revolutionary C-Scan option for MacroPATH completely automatizes this step, thus ensuring a better workflow and superior traceability.



---

AUTOMATED PROCESS

---

---

EFFORTLESS OPERATION

---

---

ENHANCED TRACEABILITY

---

---

SUPERIOR WORKFLOW

---

---


DESIGNED FOR MACROPATH

---

# DOCUMENTATION OF BIOPSIES IN A SECOND

## | C-SCAN WORKFLOW

The C-Scan option for MacroPATH enables effortless operation as it is always ready to use. The direct and error-free association between the cassette's code and its case ID offers superior traceability and therefore ultimately higher patient safety - whilst providing the pathologist with crucial information for the diagnosis.



*When receiving biopsies at grossing, the operator only needs to place the cassette under the laser beam.*



*C-Scan reads the cassette's code, creates the respective case ID folder and takes a picture at a preset zoom.*

## | DESIGNED FOR MACROPATH

MacroPATH is the Macro Digital Imaging System that addresses all the challenges related to the digitalization of the grossing step. Through the MacroPATH system, the operator can capture images, record videos and voice notes. The system includes a smart software with several assets for the grossing routine, such as measuring and sizing tools.



**MILESTONE**



HELPING  
PATIENTS

**MILESTONE Srl** - Via Fatebenefratelli, 1/5 - 24010 Sorisole (BG) - Italy  
Tel: +39 035 4128264 -  
[www.milestonemed.com](http://www.milestonemed.com) - email: [marketing@milestonemedsrl.com](mailto:marketing@milestonemedsrl.com)

**MILESTONE MEDICAL TECHNOLOGIES, INC.**  
6475 Technology Avenue, Suite F, Kalamazoo, MI 49009 - USA  
Tel: 269-488-4950 - Toll-free: 866-995-5300 -  
[www.milestonemed.com](http://www.milestonemed.com) - email: [info@milestonemed.com](mailto:info@milestonemed.com)



Distributed by:



[www.axlab.dk](http://www.axlab.dk)  
[axlab@axlab.dk](mailto:axlab@axlab.dk)  
+45 35 43 00 73