

Technical Data TFB 45		Cat. No.
Bath Capacity:	max. 1.5 liters	
Temperature:	+30 °C to +90 °C	
Rim Width:	45 mm	
Dimensions Bath Internal (Ø/H):	210 x 55 mm	
Dimensions Overall (Ø/H):	300 x 90 mm	
Weight:	ca. 2 kg	
Power Supply:	230 V / 50 – 60 Hz / 250 VA	01-4500-00
	115 V / 50 – 60 Hz / 250 VA	01-4510-00
Technical Data TFB 55		
Bath Capacity:	max. 1.5 liters	
Temperature:	+30 °C to +90 °C	
Rim Width:	55 mm	
Dimensions Bath Internal (Ø/H):	210 x 55 mm	
Dimensions Overall (Ø/H):	340 x 90 mm	
Weight:	ca. 2 kg	
Power Supply:	230 V / 50 – 60 Hz / 250 VA	01-5500-00
	115 V / 50 – 60 Hz / 250 VA	01-5510-00

Accessories	Purchasing Unit	Cat. No.
Accessories Kit for TFB 45 incl. lid, thermometer and matching holder	1 kit	11-4500-00
Accessories Kit for TFB 55 incl. lid, thermometer and matching holder	1 kit	11-5500-00
Thermometer for TFB 45 and TFB 55 , 0 – 100 °C	1 pc.	51-4551-00
Lid made of aluminum for TFB 45 and TFB 55	1 pc.	11-4540-00
Thermometer Holder for TFB 45	1 pc.	11-4541-00
Thermometer Holder for TFB 55	1 pc.	11-5541-00

MEDITE

MEDITE Medical GmbH

Wollenweberstr. 12 • 31303 Burgdorf • Germany

Phone +49 5136 88 84-0 • Fax +49 5136 88 84-55

info@medite.de • www.medite.de

Subject to change without prior notice. Pictures and colours may vary.