

Vistapath is building the next generation of pathology labs



Sentinel at a glance

Sentinel's small size and portability makes it easy to fit into any lab set up. Once implementation is complete, training lab staff is quick and simple. And best of all, it's affordable—rather than needing to find budget to buy a new tool, labs pay per sample and Vistapath manages all upkeep of the platform.



Dimensions

Sentinel module: 17" x 12" x 15" (WxDxH)

Touchscreen module: 15" x 8" x 11" (WxDxH)

Weight

Sentinel module: 15.3 lbs.

Touchscreen module: 14 lbs. including AC power brick

Operating conditions

59–82° Fahrenheit

Security

ISO 27001 certified



Electrical

Sentinel requires one 120V AC outlet

Cloud connectivity

LAN Ethernet port or LAN Wi-Fi connection

LIS integration*

Standards supported: HL7v2, CDA, X12, FHIR, Vendor APIs

Communication methods supported:

Web services, MLLP, SFTP, XDR, XDS b

*Integration supported by Redox (redoxengine.com)

Visit VISTAPATH.AI to discover how Vistapath can modernize your lab.

Created by histotechnicians and AI experts, Sentinel makes grossing fast, efficient, and easy



Vistapath's flagship product, Sentinel, is a first-of-its-kind pathology processing platform. Sentinel optimizes the grossing process using an AI-powered augmentation system that seamlessly integrates with your lab's existing workflow—helping technicians work faster and safer.

Increased accuracy and efficiency

Vistapath has taken the most time-consuming and meticulous parts of grossing and automated them—allowing technicians to complete the grossing process much faster at a higher level of quality. Sentinel photographs samples and generates a gross report of the tissue in less than one second, creating reports 90% faster. As a result, labs can process more samples in less time to increase revenue and clinical trials can quickly and effectively gather and analyze data.



Integration with current workflows

Sentinel seamlessly integrates with labs' current workflows. Technicians are directed through each grossing process, so training is quick and easy. And labs can customize parts of the workflow to make change management easy, such as the grossing mode or report templates.

About Vistapath

Vistapath was founded to accomplish a single mission: to modernize pathology by giving labs the tools they need to increase quality, consistency and efficiency in grossing workflows.



Distributed by:
Axlab A/S
+45 35 43 18 81
axlab@axlab.dk
www.axlab.dk

Axlab

YOUR VISION.
OUR INSPIRATION

Visit **VISTAPATH.AI** to discover how Vistapath can modernize your lab.