



205620 Cooling Plate CPB-2D



Distributed by:
Axlab A/S
+45 35 43 18 81
axlab@axlab.dk
www.axlab.dk

Axlab

YOUR VISION.
OUR INSPIRATION

Kunz Mobile Cooling Plate CPB-2D

Advantages

- Large cooling area
- Wheels with locking function
- Easy to move between workplaces
- Two persons can share
- Saves table space; store under the work table, lower noise
- Continuous adjustable temperature, down to -15°C
- Special height possible
- Kunz well-known product quality

Specifications

Mains power:	220-240VAC 50-60Hz
Power consumption:	250W
Temperature range at room temperature +20°C:	adjustable from +5°C to -15°C
Dimensions (WxDxH):	470x450x600mm (standard height)
Dimensions work area (WxD):	400x380mm
Weight:	23kg

