



205300 Histo Bath HIR-3D

## Kunz Histo Bath HIR-3D

### Advantages

- LED-light gives excellent visibility in the water
- Large round warm area, drying up to 27 slides
- No need for extra heating plate
- Easy to use digital touch control
- Light weight
- Easy and convenient to carry and empty the water
- An optional lid is available to cover the water
- Kunz well-known product quality

### Specifications

|                      |                             |
|----------------------|-----------------------------|
| Mains power:         | 220-240VAC 50-60Hz          |
| Power consumption:   | 190W                        |
| Temperature range:   | +20 to +80°C                |
| Temp. accuracy:      | ± 1°C                       |
| Internal dimensions: | Diameter 217mm, Depth 61mm  |
| External dimensions: | Diameter 323mm, Height 88mm |
| Weight:              | 1,5kg                       |



[www.kunzinstruments.se](http://www.kunzinstruments.se)   [info@kunzinstruments.se](mailto:info@kunzinstruments.se)



Distributed by:  
**Axlab A/S**  
+45 35 43 18 81  
[axlab@axlab.dk](mailto:axlab@axlab.dk)  
[www.axlab.dk](http://www.axlab.dk)

**Axlab**

YOUR VISION.  
OUR INSPIRATION